

Product Datasheet

GSTA1 antibody (orb167005)

Description	Mouse monoclonal antibody to GSTA1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	GSTA1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM,k
Clonality	Monoclonal
MW	25631
Uniprot ID	P08263
NCBI	NP_665683.1
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.

