

Product Datasheet

ATG4C antibody (orb167006)

Description	Mouse monoclonal antibody to ATG4C.
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ATG4C
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	52497
Uniprot ID	Q96DT6
NCBI	NP_835739.1 , NP_116241.2
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.

