

Product Datasheet

CD11c antibody (orb167010)

Description	Mouse monoclonal antibody to CD11c
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ITGAX
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	127829
Uniprot ID	P20702
NCBI	NP_000878.2
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

