

Product Datasheet

AKT1S1 antibody (orb167011)

Description	Mouse monoclonal antibody to AKT1S1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	AKT1S1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	27383
Uniprot ID	Q96B36
NCBI	NP_001092103.1 , NP_001265089.1 , NP_001092102.1 , NP_115751.3 , NP_001265088.1
Dilution Range	WB: 1:100
Expiration Date	12 months from date of receipt.

