

Product Datasheet

SWAP70 antibody (orb167015)

Description	Mouse monoclonal antibody to SWAP70.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	SWAP70
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	68998
Uniprot ID	Q9UH65
NCBI	NP_055870.2
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

