




Product Datasheet

PSMA5 antibody (orb167018)

Description	Mouse monoclonal antibody to PSMA5.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	PSMA5
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	26411
Uniprot ID	P28066
NCBI	NP_002781.2 , NP_001186701.1 , NP_001186702.1 , NP_001186703.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.



WB analysis of (lane 1) human liver lysa...



Immunofluorescence analysis of Hela cell...