

Product Datasheet

Goat Estradiol (E2) ELISA Kit (orb1670188)

Description	Goat Estradiol (E2) ELISA Kit
Reactivity	Goat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 40-1400 pg/mL. Sensitivity: 18 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.