

Product Datasheet

RUVBL1 antibody (orb167019)

Description	Mouse monoclonal antibody to RUVBL1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RUVBL1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50228
Uniprot ID	Q9Y265
NCBI	NP_003698.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

