

Product Datasheet

C-reactive protein antibody (orb167021)

Description	Mouse monoclonal antibody to C-reactive protein
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CRP
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
MW	25039
Uniprot ID	P02741
NCBI	NP_000558.2
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

