

Product Datasheet

CD57 antibody (orb167022)

Description	Mouse monoclonal antibody to CD57
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	B3GAT1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	38256
Uniprot ID	Q9P2W7
NCBI	NP_061114.2 , NP_473366.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

