



Product Datasheet

GAD2 antibody (orb167024)

Description	Mouse monoclonal antibody to GAD2.
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GAD2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	65411
Uniprot ID	Q05329
NCBI	NP_000809.1 , NP_001127838.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...