

Product Datasheet

Pig Glial fibrillary acidic protein (GFAP) ELISA Kit (orb1670290)

Description	Pig Glial fibrillary acidic protein (GFAP) ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.31-20 ng/mL. Sensitivity: 0.156 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.