

## Product Datasheet

### MX1 antibody (orb167030)

|                            |                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to MX1                                                                                                                                    |
| <b>Species/Host</b>        | Mouse                                                                                                                                                               |
| <b>Reactivity</b>          | Human                                                                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                        |
| <b>Tested Applications</b> | WB                                                                                                                                                                  |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                   |
| <b>Target</b>              | MX1                                                                                                                                                                 |
| <b>Form/Appearance</b>     | Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS. |
| <b>Note</b>                | For research use only                                                                                                                                               |
| <b>Isotype</b>             | IgG2b                                                                                                                                                               |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                          |
| <b>MW</b>                  | 75520                                                                                                                                                               |
| <b>Uniprot ID</b>          | <a href="#">P20591</a>                                                                                                                                              |
| <b>NCBI</b>                | <a href="#">NP_001171517.1</a> , <a href="#">NP_001138397.1</a> , <a href="#">NP_002453.2</a>                                                                       |
| <b>Dilution Range</b>      | WB: 1:100-500                                                                                                                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                     |

