

Product Datasheet

Human Probable G-protein coupled receptor 116 (GPR116) ELISA Kit (orb1670307)

Description Human Probable G-protein coupled receptor 116 (GPR116) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.156-10 ng/mL. Sensitivity: 0.078 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.