

Product Datasheet

CD13 antibody (orb167031)

Description	Mouse monoclonal antibody to CD13
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ANPEP
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	109540
Uniprot ID	P15144
NCBI	NP_001141.2
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

