

Product Datasheet

Fish Heat Shock Protein 90 (HSP-90) ELISA Kit (orb1670360)

Description	Fish Heat Shock Protein 90 (HSP-90) ELISA Kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.5-8 ng/mL. Sensitivity: 0.255 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.