

Product Datasheet

Human Hydroxyproline (Hyp) ELISA Kit (orb1670368)

Description Human Hydroxyproline (Hyp) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 62.5-4000 ng/mL.
Sensitivity: 35.5 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.