

Product Datasheet

Human Interleukin-10 receptor subunit beta (IL10RB) ELISA Kit (orb1670414)

Description	Human Interleukin-10 receptor subunit beta (IL10RB) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 31 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.