

Product Datasheet

Rabbit Interleukin 2 (IL-2) ELISA Kit (orb1670437)

Description	Rabbit Interleukin 2 (IL-2) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 16.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.