

Product Datasheet

MC1R antibody (orb167046)

Description	Mouse monoclonal antibody to MC1R
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	MC1R
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	34706
Uniprot ID	Q01726
NCBI	NP_002377.4
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

