

Product Datasheet

OR5B12 antibody (orb167049)

Description	Mouse monoclonal antibody to OR5B12 (C-term).
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	OR5B12
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Ig2b
Clonality	Monoclonal
MW	35182
Uniprot ID	Q96R08
NCBI	NP_001004733.1
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

