




Product Datasheet

ITGB2 antibody (orb167054)

Description	Mouse monoclonal antibody to ITGB2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ITGB2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	84782
Uniprot ID	P05107
NCBI	NP_000202.2 , NP_001120963.1
Dilution Range	IHC-P: 1:25, WB: 1:100-500
Expiration Date	12 months from date of receipt.



WB analysis of HL-60 cell line using ITG...



Immunohistochemical analysis of formalin...