

Product Datasheet

IGHD antibody (orb167055)

Description	Mouse monoclonal antibody to IGHD (N-term).
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	IGHD
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Euglobin precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
MW	42353
Uniprot ID	P01880
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

