

Product Datasheet

Human Kynurenine (KYN) ELISA Kit (orb1670559)

Description	Human Kynurenine (KYN) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.12-500 pmol/mL. Sensitivity: 0.12 pmol/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.