

Product Datasheet

Human Laminin subunit alpha-2 (LAMA2) ELISA Kit (orb1670567)

Description Human Laminin subunit alpha-2 (LAMA2) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 3.12-200 ng/mL. Sensitivity: 1.34 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.