

Product Datasheet

DRD4 antibody (orb167061)

Description	Mouse monoclonal antibody to DRD4
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DRD4
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43901
Uniprot ID	P21917
NCBI	NP_000788.2
Dilution Range	WB: 1:100-250
Expiration Date	12 months from date of receipt.

