

Product Datasheet

Fish Leptin (LEP) ELISA Kit (orb1670611)

Description	Fish Leptin (LEP) ELISA Kit
Reactivity	Fish
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 62.5-1000 pg/mL. Sensitivity: 31.25 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.