

Product Datasheet

Angiotensinogen antibody (orb167062)

Description	Mouse monoclonal antibody to Angiotensinogen
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AGT
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53154
Uniprot ID	P01019
NCBI	NP_000020.1
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

