

Product Datasheet

CYP1A1 antibody (orb167063)

Description	Mouse monoclonal antibody to CYP1A1 (C-term).
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CYP1A1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	58165
Uniprot ID	P04798
NCBI	NP_000490.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

