




Product Datasheet

PINK1 antibody (orb167066)

Description	Mouse monoclonal antibody to PINK1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	PINK1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	62769
Uniprot ID	Q9BXM7
NCBI	NP_115785.1
Dilution Range	WB: 1:100-500, IHC-P: 1:50-100, IF/ICC: 1:25
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Immunohistochemical analysis of formalin...