

Product Datasheet

Mouse Lysophosphatidylcholine acyltransferase 1 (LPCAT1) ELISA Kit (orb1670677)

Description	Mouse Lysophosphatidylcholine acyltransferase 1 (LPCAT1) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.056 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.