

Product Datasheet

Mouse Tumor necrosis factor receptor superfamily member 3 (LTBR) ELISA Kit (orb1670727)

Description Mouse Tumor necrosis factor receptor superfamily member 3 (LTBR) ELISA Kit

Reactivity Mouse

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 31.25-2000 pg/mL. Sensitivity: 7.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.