




Product Datasheet

EGFR antibody (orb167074)

Description	Mouse monoclonal antibody to EGFR.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	EGFR
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1κ
Clonality	Monoclonal
MW	134277
Uniprot ID	P00533
NCBI	NP_958439.1 , NP_005219.2 , NP_958441.1 , NP_958440.1
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of HepG2 cells ...