

Product Datasheet

Human IFN- α ELISA KIT (orb1670786)

Description	Human IFN- α ELISA KIT
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.625-1000 pg/mL. Sensitivity: 7 pg/mL. Assay Time: 3 hours. Sample Volume: 20 μ L. Sample Type: Cell culture supernatant, serum, plasma (EDTA, citrate, heparin). Detection method: Colorimetric.
Expiration Date	Please enquire.