

Product Datasheet

IKK gamma antibody (orb167080)

Description	Mouse monoclonal antibody to IKK gamma.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	IKBKG
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	48198
Uniprot ID	Q9Y6K9
NCBI	NP_001093326.2 , NP_003630.1 , NP_001138727.1 , NP_001093327.1
Dilution Range	WB: 1:100-500
Expiration Date	12 months from date of receipt.

