




Product Datasheet

GFP Tag antibody (orb167086)

Description	Mouse monoclonal antibody to GFP Tag.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	Purified His-tagged GFP protein was used to produced this monoclonal antibody.
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	27 k
Uniprot ID	C5MKY7
Dilution Range	WB: 1:4000, IF/ICC: 1:50-100, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of GFP recombinant...



Confocal immunofluorescence analysis of ...