

Product Datasheet

Rabbit apolipoprotein B100 (Apo-B100) ELISA Kit (orb1670865)

Description	Rabbit apolipoprotein B100 (Apo-B100) ELISA Kit
Reactivity	Rabbit
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78 ng/mL-5000 ng/mL. Sensitivity: 19.5 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.