

Product Datasheet

Plant Vitamin A,VA ELISA Kit (orb1670867)

Description	Plant Vitamin A,VA ELISA Kit
Reactivity	Plant
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625 ng/mL-10 ng/mL. Sensitivity: 0.467 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: plant tissue.
Expiration Date	Please enquire.