

Product Datasheet

Pig progesterone,PROG ELISA Kit (orb1670872)

Description	Pig progesterone,PROG ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.2 ng/mL-80 ng/mL. Sensitivity: 0.2 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.