

Product Datasheet

CD3 antibody (orb167090)

Description	Mouse monoclonal antibody to CD3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CD3E
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23147
Uniprot ID	P07766
NCBI	NP_000724.1
Dilution Range	FACS: 1:100
Expiration Date	12 months from date of receipt.

