

Product Datasheet

Monkey Progesterone(PROG) ELISA Kit (orb1670933)

Description	Monkey Progesterone(PROG) ELISA Kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.25 ng/mL-100 ng/mL. Sensitivity: 0.25 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 μ L. Sample Type: serum, plasma, cell culture supernates, tissue homogenates.
Expiration Date	Please enquire.