

Product Datasheet

Monkey luteinizing hormone(LH)ELISA kit (orb1670934)

Description	Monkey luteinizing hormone(LH)ELISA kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 2 mIU/mL-75 mIU/mL. Sensitivity: 0.5 mIU/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.