

Product Datasheet

Monkey estradiol (E2) ELISA kit (orb1670936)

Description	Monkey estradiol (E2) ELISA kit
Reactivity	Monkey
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 40 pg/mL-1000 pg/mL. Sensitivity: 25 pg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.