

## Product Datasheet

# Human kappa immunoglobulin light chain,κ-IgLC ELISA Kit (orb1670987)

|                            |                                                                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human kappa immunoglobulin light chain,κ-IgLC ELISA Kit                                                                                                              |
| <b>Reactivity</b>          | Human                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                |
| <b>Application notes</b>   | Detect Range: 0.078 µg/mL-20 µg/mL. Sensitivity: 0.046 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                      |