

## Product Datasheet

# Human anti-Chlamydia trachomatis(CT) antibody (IgG) ELISA Kit (orb1671065)

|                            |                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human anti-Chlamydia trachomatis(CT) antibody (IgG) ELISA Kit                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA                                                                                                                                                               |
| <b>Note</b>                | For research use only                                                                                                                                               |
| <b>Application notes</b>   | Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                     |