



Product Datasheet

DRD3 antibody (orb167107)


Description	Rabbit polyclonal antibody to DRD3.
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Dopamine Receptor D3 (Dopamine Receptor D3 (DRD3))(NP_000787.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	P35462
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of Mou...



Immunohistochemical staining of rat brai...



Immunofluorescence analysis of U20S cell...