

Product Datasheet

Horse tri-iodothyronine (T3) ELISA kit (orb1671079)

Description	Horse tri-iodothyronine (T3) ELISA kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.5 ng/mL-8 ng/mL. Sensitivity: 0.25 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma, tissue homogenates.
Expiration Date	Please enquire.