

Product Datasheet

Horse Progesterone(PROG) ELISA kit (orb1671082)

Description	Horse Progesterone(PROG) ELISA kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0 ng/mL-80 ng/mL. Sensitivity: 0.246 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.