

Product Datasheet

Horse Albumin(Alb) ELISA Kit (orb1671087)

Description	Horse Albumin(Alb) ELISA Kit
Reactivity	Equine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.25 µg/mL-80 µg/mL. Sensitivity: 1.247 µg/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.