

Product Datasheet

Guinea Pig Cortisol ELISA Kit (orb1671095)

Description	Guinea Pig Cortisol ELISA Kit
Reactivity	Guinea pig
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.049 ng/mL-200 ng/mL. Sensitivity: 0.049 ng/mL. Assay Time: 1-5 h. Detection wavelength: 450 nm. Sample Volume: 50-100 µL. Sample Type: serum, plasma.
Expiration Date	Please enquire.